



ENGLISH

BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

INTENDED USER PROFILE:

THE USER SHOULD BE ABLE TO:
• READ AND UNDERSTAND THE DIRECTIONS,
WARNINGS AND CAUTIONS.

INTENDED USE/INDICATIONS: Minor knee and patellofemoral (knee cap) instability.

CONTRAINDICATIONS: None.

WARNINGS AND PRECAUTIONS:

- If you experience any pain, swelling, sensation changes, or any unusual reactions while using this product, consult your medical professional immediately.
- **The Patella Donut contains Neoprene.** Neoprene supports should not be worn by individuals with known susceptibility to dermatitis or allergies to neoprene by-products.

APPLICATION INFORMATION:

1. Open and detach anterior (front) and posterior (back) straps. In seated position, pull brace on.
2. Centre brace over knee cap.
3. Pull each anterior (front) strap through D-Ring.
4. Secure each strap on the velcro.
5. Repeat steps 3 and 4 for posterior (back) straps.
6. Adjust anterior (front) straps to appropriate comfort level.
7. Repeat step 6 for posterior (back) straps.
8. To trim donut, remove hinges, turn brace inside out.
 - Determine support / shape required (Superior, Inferior, Lateral or Medial)
 - Find best opening in donut cover
 - Cut the donut with scissors
 - Remove excess donut material not required

USE AND CARE: Hand wash in water (30°) using mild soap. Air dry only.

WARRANTY: DJO, LLC will repair or replace all or part of the unit and its accessories for material or workmanship defects for a period of six months from the date of sale.

INTENDED FOR SINGLE PATIENT USE
NOT MADE WITH NATURAL RUBBER LATEX

NOTICE: WHILE EVERY EFFORT HAS BEEN MADE IN STATE-OF-THE-ART TECHNIQUES TO OBTAIN THE MAXIMUM COMPATIBILITY OF FUNCTION, STRENGTH, DURABILITY AND COMFORT, THERE IS NO GUARANTEE THAT INJURY WILL BE PREVENTED THROUGH THE USE OF THIS PRODUCT.

ESP

ANTE:
TES IN
LA CO
IMPRE
PERFIL
EL USU
- LEER Y
Y PREC

USO P
menor

CONTI

ADVEF

- Si exp
sensai
produ
- El anil
persor
alergia
neopre

INSTRU

1. Abra y
y post
cológi
 2. Centre
 3. Tire de
anilla e
 4. Fije ca
 5. Repita
posteri
 6. Ajuste
que le
 7. Repita
(trase
 8. Para rec
la rodill
- Deter
(super
 - Locali
anillo e
 - Corte e
 - Quite e

USO Y CUI
con un jab

GARANTÍA
unidad o pe
material del
durante el p
de venta.

REV A

© 2016 DJO, LLC

EU Authorized Representative
MDSS GmbH
Schiffgraben 41
30175 Hannover
Germany
Phone (+49) 511-6262 8630
Phone (+49) 511-6262 8633

EC REP

DJO Global Company
1430 Decision Street
Vista, CA 92081-8553 • USA
T 760.727.1280 • F 760.734.5608
DJOglobal.com



Quick Fit Hinged Knee Brace

RODILLERA ARTICULADA QUICK FIT

QUICK FIT KNEIORTHESE MIT GELENK

GINOCCHIERA SNODATA QUICK FIT

ATTELLE DE GENOU ARTICULÉE QUICK FIT

QUICK FIT KNEIBRACE MET SCHARNIER

DONJOY

INDICADO PA
NO FABRICAD

AVISO: SI BIEN SE HA
TÉCNICAS MÁS MOD
DE LA FUNCIÓN, REI
GARANTÍA ALGUNAI
DE ESTE PRODUCTO.

