
Esteem™ Soft ConvaTec

Instructions for use /Bruksanvisning / Gebrauchsanleitung / Návod k použití / *Návod na použitie / Instruccia stosowania / Instruções por uso / OAHŤIEĽ XPŤHZE / Brugsanvisning / Käyttöohjeet / Bruksanvisning / Kullannin Talimat / Naudojimo instrukcija / Uputstvo za upotrebu / Upute za uporabu / Lietošanas pamācība / Navodila za uporabo / Kasutusjuhend / Használati útmutató / Gebrauchsanleitung / Notice d'utilisation / طريقة الإستعمال / Instrucciones de uso / Instruções de utilização / Instruções de utilização*

ConvaTec

ENGLISH

Please read these instructions for use carefully, as they contain important information about the use of our Esteem™ + Soft ConvaTec products.

PRODUCT DESCRIPTION
Esteem™+ Soft ConvaTec pouches are one-piece closed, drainable or urostomy pouches with a pre-cut or cut-to-fit skin barrier and an internal adhesive.

INTENDED USE: Management of stoma output.

INDICATION FOR USE: For use with a colostomy, ileostomy or urostomy.

PRECAUTIONS AND OBSERVATIONS
This device is for single-use only and should not be re-used. Please may lead to increased risk of infection or cross contamination. Physical properties of the device may never be optimal for intended use.

If during the use of this device or as a result of its use a serious incident has occurred, please report it to the manufacturer and to your national authority. For cut-to-fit pouches, make sure the opening you cut is not too large. An ill-fitting pouch can cause skin irritation and inflammation. Be very careful not to increase the opening size past the outer ring marked on the water or to the edge of or exposing the inner plastic convex disc.

Note that different pouch types are appropriate for different stomas. Always consult your healthcare professional on the most appropriate pouch type for your stoma prior to use.

Important: A healthcare professional should assess the need for skin barriers with convexity prior to use.
DIRECTIONS FOR USE
1. Preparation
Clean the skin around your stoma with water and a soft cloth. Pat the area dry. Do not use oil-based products as they may interfere with the adhesion of the skin barrier.

Please make sure you know the diameter of your stoma. Esteem™ + Soft ConvaTec is available in pre-cut and cut-to-fit options.

For cut-to-fit pouches you can use the cutting guide on the water hereafter provided. The hole in the adhesive water should fit the stoma as closely as possible. Cut the stoma hole in the water to the diameter of your stoma or a little bit smaller (figure 1).
2. Use of the Esteem™ + Soft ConvaTec system
Partially inflate the pouch by blowing to separate the inner layers of the pouch. Warm the skin barrier with the palm of your hand for a few seconds to aid adhesion. Remove the plastic covering from the skin barrier (figure 2). Try to avoid touching the adhesive area too much with your hands.

Hold the skin barrier below your stoma. Beginning from the bottom carefully adhere the barrier to the skin around your stoma (figure 3).
Smooth over the back of the adhesive area with your fingers and press gently but firmly around your stoma (figure 4).

3. Advice for use of a drainable pouch with Safe Seal™ closure
For those with an ileostomy or colostomy there is a type of pouch that can be drained, the tail closure cleaned and sealed again.

Please note the following advice for using a drainable pouch with Safe Seal™ closure:

Make sure the bottom of the pouch is securely fastened. Fold up the neck of the pouch and continue to use for a total of 5 folds until the two parts of the fastener are in contact with each other. Press firmly together to gain a secure seal (figure 5).

Apply the pouch as instructed. Pull the tab to open the pouch and empty the pouch.

4. Advice for using a urostomy pouch
Make sure the bottom of the pouch is securely fastened. Apply the pouch as instructed. Use the supplied adaptor to connect the urostomy pouch to a night drainage pouch or leg bag if necessary (figure 6).

5. Information about pouches with filters
There is a built in filter in closed and drainable pouches through which gases are absorbed and neutralized. The filter either provides a passive or active protection against moisture (e.g. when bathing, showering) and ensures thereby the functional performance of the filter.

Cover the filter with the sticker when bathing, swimming, etc (figure 7). Remove after water activity.

6. Pouch removal and disposal
To remove an old pouch for disposal, gently peel downwards while holding the skin above the water with your free hand (figure 8). This helps avoid skin irritation and makes the removal easier.

Empty the contents of the pouch into the toilet.
Never flush the pouch down the toilet.
After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local laws and regulations.

Made in Germany
© 2021 ConvaTec Inc.
All trademarks are the property of their respective owners. ESTEEM is a registered trademark of ConvaTec Inc. in the United States.

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ARTWORK SPECIFICATION	
Artcode	SI2358VA
Artwork Description	INS EST+ SOFT CVX EUR
ICC	422527, 422528, 422529, 422530, 422531, 422532, 422533, 422534, 422535, 422536, 422537, 422538, 422539, 422540, 422541, 422542, 422543, 422544, 422545, 422546, 422547, 422548, 422549, 422550, 422551, 422552, 422553, 422554, 422555, 422556, 422557, 422558, 422559, 422560, 422561, 422562, 422563, 422564, 422565, 422360, 422361, 422362, 422363, 422364, 422365
SAP	1727310, 1727312, 1727314, 1727316, 1727318, 1727320, 1727322, 1727324, 1727326, 1727328, 1727330, 1727332, 1727334, 1727336, 1727338, 1727340, 1727342, 1727344, 1727346, 1727348, 1727350, 1727352, 1727354, 1727356, 1727358, 1727360, 1727362, 1727364, 1727366, 1727368, 1727370, 1727372, 1727374, 1727376, 1727378, 1727380, 1727382, 1727384, 1727386, 1724312, 1724313, 1724314, 1724318, 1724319, 1724320
KOM/REF	N/A
IS Code	N/A
Component Type	Insert
LR#	LR-004354
Existing Art Code being Revised	SI2358
REQUIREMENT	SPECIFICATION
Bucket	EUR
Local Code(s)	N/A
Size Code / Dieline	N/A
Dimensions	Flat: 24'' x 28'' (609.6mm x 711.2mm) Folded: 4'' x 5.6'' (101.6mm x 142.2mm)
Color(s) / Colour(s)	Black
Specification Reference(s)	TF23-048, TP23-048
Bar Code: UPC (A) or GTIN*	N/A
Bar Code: PZN – Code 39*	N/A
Bar Code: JAN-13*	N/A
* NOTE: Bar Code: See artwork for positioning, magnification and truncation. Scanability verified to ANSI grade “C” or better.	
HISTORY	
6-June-2018 / SI2358: New insert for Esteem+ Soft Convex. 3-Apr-2019 / SI2358: IFU utilization extended for new ICC and finished codes. 9-JUNE-2021 / SI2358VA: All symbols updated in accordance with ISO 15223-1:2016 and MDR compliance. Moved caution symbol together with the rest of the symbols. Added language codes. Added MD symbol and the corresponding definition in all languages required. Added statements in the IFU Intended use and indication for use and it's translation acc. to SM-076885. Updated IFU statements regarding the device single use, during and after use according to SM-076885. Art version changed to VA. Copyright and revision date updated. Convatec logo updated. Updated Deeside address and manufacturing symbol. Added importers symbol and explanation. EC REP symbol with address added with importers symbol.	