


PROVOX[®] LaryTube[™]

MD

**Rx
ONLY**

CE 0413

42° C
108° F
RT
2° C
36° F

Figure 1


Figure 2


Figure 3


Figure 4


Figure 5


Figure 6a


Figure 6b


Figure 6c


Figure 6d


Figure 6e


Figure 6f


Figure 6g



Figure 6h

Prescription information

CAUTION: United States Federal law restricts this device to sale, distribution and use by or on order of a physician or a licensed practitioner. The availability of this product without prescription outside the United States may vary from country to country.

Disclaimer

Atos Medical offers no warranty - neither expressed nor implied - to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

Patents and trademarks

Provox[®] is a registered trademark owned by Atos Medical AB, Sweden. LaryTube[™] and LaryClip[™] are trademarks of Atos Medical AB.

For information about protective rights (e.g. patents), please refer to our web page www.atosmedical.com/patents.

Manufacturer; Hersteller; Fabrikant; Fabricant; Produttore;
 Fabricante; Fabricante; Tillverkare; Producent; Produsent;
 Valmistaja; Framleiðandi; Tootja; Ražotājs; Gamintojas;
 Výrobce; Gyártó; Výrobca; Proizvajalec; Producent; Producător;
 Proizvođač; Proizvođač; Κατασκευαστής; Производител; Üretici;
 მწარმოებელი; Ἀρτοποιητής; İstehsalçı; Производитель; Pro-
 dusen; Pengilang; 製造業者; 제조사; 製造商; 制造商;
 היצרן; الجهة المصنعة

Date of manufacture; Herstellungsdatum; Productiedatum; Date de
 fabrication; Data di fabbricazione; Fecha de fabricación; Data de
 fabrico; Tillverkningsdatum; Fremstillingsdato; Produksjonsdato
 ; Valmistuspäivä; Framleiðsludagsetning; Valmistamise kuupäev;
 Ražošanas datums; Pagaminimo data; Datum výroby; Gyártás
 időpontja; Dátum výroby; Datum proizvodnje; Data produkcji; Data
 fabricației; Datum proizvodnje; Datum proizvodnje; Ημερομηνία
 κατασκευής; Дата на производство; Üretim tarihi; წარმოების
 თარიღი; Ἀρτοποιημάτων αμύση; İstehsal tarixi; Дата
 изготовления; Tanggal produksi; Tarikh pembuatan; 製造日;
 제조일자; 製造日期; 制造日期; تاريخ التصنيع

Use-by date; Verfallsdatum; Houdbaarheidsdatum; Date limite
 d'utilisation; Data di scadenza; Fecha de caducidad; Data de vali-
 dade; Används senast; Anvendes inden; Brukes innen; Viimeinen
 käyttöpäivä; Síðasti notkunardagur; Kölblik kuni; Derīguma
 termiņš; Naudoti iki; Datum expirace; Felhasználatoság dátuma;
 Dátum najneskoršej spotreby; Rok uporabnosti; Data ważności; Data
 expirării; Rok uporabe; Rok upotrebe; Ημερομηνία λήξης; Срок на
 годност; Son kullanım tarihi; ვარგისობა; Ἀχρησιμοποίητων
 δαμύλη; Son istifadə tarixi; Срок годности; Gunakan sebelum

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
1. Descriptive information

1.1 Intended use

The Provox LaryTube is a holder for devices in the Provox HME System intended for vocal and pulmonary rehabilitation after total laryngectomy.

For patients with a shrinking tracheostoma it is also used to maintain the tracheostoma for breathing.

The Provox LaryTube is intended for single patient use.

1.2 CONTRAINDICATIONS

Provox LaryTube is not intended to be used by patients that:

- are under any form of mechanical ventilation.
- have damaged tracheal or tracheostoma tissue.

1.3 Description of the device

The Provox LaryTube is a tube made of medical grade silicone rubber. The purpose of the device is to create a comfortable and airtight fit between the Provox LaryTube and the tracheostoma, and also to provide attachment for devices from the Provox HME System.

They are delivered single packed, non-sterile.

There are 3 models: Standard, Fenestrated and with Ring (Fig. 1).

The Standard model and the Ring version can be fenestrated so that air can go through the voice prosthesis for voice prosthesis users. The holes are punched by using the Provox Fenestration Punch according to the Instructions for Use accompanying the Provox Fenestration Punch, see ordering information.

Standard versions (Fig. 1a) – made for use with or without a voice prosthesis. Can be attached with a Provox TubeHolder or Provox LaryClips.

Fenestrated versions (Fig. 1b) – For voice prosthesis users. Can be attached with a Provox TubeHolder or Provox LaryClips.

Ring versions (Fig. 1c) – made for use with or without a voice prosthesis. Can only be attached with a Provox Adhesive.

The different parts of the Provox LaryTube are (Fig. 1d-i):

- d) Tube
- e) Shield (conical)

- f) HME and Accessory Holder
- g) Wings
- h) Ring for adhesive connection
- i) Fenestration holes

1.4 WARNINGS

- **DO NOT** reuse a used Provox LaryTube in another patient. The device is intended for single patient use only. Reuse in another patient may cause cross contamination.
- **DO NOT** use lubricating gels if the patient has allergies related to these substances.
- **DO** ensure that the patient has been trained in the use of the device. The patient should demonstrate the ability to understand and consistently follow the Instructions for Use without clinician supervision.

1.5 PRECAUTIONS

- Always use a Provox LaryTube of an appropriate size. If the Provox LaryTube is too wide, too narrow, too long or too short it may cause tissue damage, bleeding or irritation. In addition, it may be harder to produce voice, because the LaryTube may obstruct the voice prosthesis. Inserting, removing and pressing on the HME Cassette may exert pressure on the voice prosthesis if the Provox LaryTube is not sized appropriately. A LaryTube that is too narrow may cause the stoma to shrink.
- **DO NOT** use force during insertion. Always insert and remove the Provox LaryTube according to the instructions below. If the product is pushed too far into trachea it may obstruct breathing. The use of force during insertion may cause tissue damage, bleeding or irritation and accidental dislodgement of the voice prosthesis and/or XtraFlange if used.
- If fenestration is necessary, make sure that it is done with the Provox Fenestration Punch. This will ensure that the holes are small and smooth. Fenestration holes that are too large may cause dislodgement of the voice prosthesis or XtraFlange. Fenestration holes with sharp edges may cause formation of granulation tissue. Make sure that no silicone residues are left in the holes or lumen of the LaryTube after fenestration.
- **DO NOT** use a dirty or contaminated Provox LaryTube. Clean and disinfect the device according to the cleaning and disinfection instructions below.
- **DO NOT** use a damaged device since this might cause aspiration of small parts or cause formation of granulation tissue.
- Only use water-soluble lubricants. Oil-based lubricants (e.g., Vaseline) may weaken, damage, or destroy the product and should not be used.
- Carefully monitor the tracheostoma tissue during radiation therapy. Stop using the device if the stoma gets irritated or starts to bleed.

- Patients with bleeding disorders or patients on anticoagulants should not use the device if it causes bleeding that re-occurs.
- Stop using the device if granulation tissue forms in the tracheostoma.
- **DO NOT** clean the device while inside the stoma since this could cause tissue damage. Always remove the device from the stoma before cleaning.

2. Instructions for use

2.1 Preparation

Choosing correct size

Diameter: The Provox LaryTubes are available in 4 different diameters, 8, 9, 10 and 12. Measure the size of the tracheostoma (largest diameter) with a ruler and fit the Provox LaryTube accordingly. In irregularly shaped tracheostomas, it may be difficult to achieve an airtight fit.

Length: The Provox LaryTubes are available in 3 different lengths, 27, 36 and 55 mm.

To help choose the size(s) of LaryTube to prescribe to the patient, the clinician can use a Provox LaryTube Sizer Kit. The Sizer Kit contains samples (“sizers”) of commercially available Provox LaryTubes.

CAUTION: During fitting always make sure that the LaryTube does not catch on to the tracheal flange of the voice prosthesis (Fig. 3).

Preparing for insertion

Before first use check that you have the correct size and that the package is not damaged or opened. Do not use otherwise.

Make sure that your hands are thoroughly clean and dry before placing them in the area of the stoma or handling the Provox LaryTube and/or accessories.

Have a mirror and sufficient light to illuminate the stoma.

Always carefully inspect the Provox LaryTube **BEFORE EACH USE** (i.e. before insertion), and ensure that it is undamaged and shows no signs of tears, cracks, or crusts. If you detect any damage, do not use the product, and obtain a replacement.

Before insertion, rinse the device with drinking water to remove any disinfectant residues (Fig. 2). Make sure that you only insert clean devices. In a hospital setting it is better to use sterile water rather than drinking water.

If needed, lightly lubricate the Tube with a small amount of water-soluble lubricant. The following water-soluble lubricants are recommended for use with the Provox LaryTube:

- KY Jelly[®]
- SurgiLube[®]
- Xylocaine[®]/(lidocaine)

CAUTION: Only use water-soluble lubricants. Oil-based lubricants (e.g.,

Vaseline) may weaken, damage, or destroy the product and should not be used at any time.

Make sure that you do not lubricate the HME Holder, HME Cassette or any accessory which is to be held by the Provox LaryTube. This may lead to accidental detachment.

2.2 Operating instruction

Insertion (Fig. 3):

Provox LaryTube Standard and Fenestrated:

Gently insert the Provox LaryTube directly into the tracheostoma. Attach these tubes with either Provox LaryClip or Provox TubeHolder (Fig. 4b-c).

Provox LaryTube with Ring:

Attach a Provox Adhesive and then gently insert the Provox LaryTube through the adhesive into the tracheostoma (Fig. 4a).

Make sure that the adhesive is centered properly so that the LaryTube does not press on the sides of the stoma.

Insert the Provox HME System component into the HME Holder.

NOTE: Some patients might experience more coughing after inserting the LaryTube for the first time or when occluding the HME to speak. This usually diminishes as the user adjusts to the device. If discomfort persists the LaryTube should be removed.

Device removal

1. Hold the Provox LaryTube in place and remove the Provox HME System component from the tube (Fig. 5).
2. Remove the Provox LaryTube from the stoma.

NOTE: Insertion and/or removal of the LaryTube may occasionally cause some slight bleeding, irritation or coughing. Inform your clinician if these symptoms persist.

2.3 Cleaning and disinfection

CAUTION: DO NOT clean the device while inside the stoma. This can cause injury. Always remove the device from the stoma before cleaning.

Clean the Provox LaryTube after each use. Do not use any water other than drinking water to clean and rinse the device.

CAUTION: Be careful not to burn your fingers in the warm cleaning water.

1. Rinse the Provox LaryTube in drinking water (Fig. 6a).

2. Place the Provox LaryTube in non-oil based dish soap (Fig. 6b).
3. Clean the inside of the tube with the Provox TubeBrush (Fig. 6c, 6h).
4. If a fenestrated tube is used the holes should be cleaned with a Provox Brush, see ordering information (Fig. 6d, 6h).
5. Rinse the Provox LaryTube in drinking water (Fig. 6e).
6. Disinfect the Provox LaryTube at least once a day (Fig. 6f) with one of the following methods:
 - Ethanol 70% for 10 minutes
 - Isopropylalcohol 70% for 10 minutes
 - Hydrogenperoxide 3% for 60 minutes

CAUTION: Do not use the device until it is completely dry (Fig. 6g). Inhalation of disinfectant fumes can cause severe coughing and irritation of the airways.

If the Provox LaryTube looks dirty or has air dried in an area with a risk of contamination the device should be both cleaned and disinfected before use. A risk of contamination could exist if the device has been dropped on the floor or if it has been in contact with a pet, someone with a respiratory infection, or any other gross contamination.

During hospitalization it is important to both clean and disinfect the Provox LaryTube after use but also before use since there is an increased risk of infection and contamination. In a hospital it is better to use sterile water for cleaning and rinsing, rather than drinking water.

CAUTION: Do not clean or disinfect by any other method than described above as this might cause product damage and patient injury.

2.4 Storage instructions

When not in use, clean and disinfect the device as described above and then store the Provox LaryTube in a clean and dry container at room temperature. Protect from direct sunlight.

2.5 Device lifetime

The Provox LaryTube may be used for a maximum period of 6 months. If the device shows any signs of damage it shall be replaced earlier.

2.6 Disposal

Always follow medical practice and national requirements regarding biohazard when disposing of a used medical device.

2.7 Accessories

CAUTION: Use only original components. Other accessories may cause product damage and/or malfunction.

Some of the devices that can be used with the Provox LaryTube are:

- Provox HME Cassettes (Heat and Moisture Exchangers)
- Provox FreeHands HME, for speaking without manual stoma occlusion.
- Provox ShowerAid that can help to prevent water from entering the LaryTube during showering.
- A Provox TubeHolder or Provox LaryClips can be used for additional support to prevent the Provox LaryTube from falling out of the stoma.
- For Provox LaryTube with Ring a Provox Adhesive is used as attachment and keeping the LaryTube in place.

Please see ordering information for different accessories available.

3. Troubleshooting information

Bleeding of the stoma

Stop using the device. After the tracheostoma is healed gradually try to increase the duration of time when using the device. If even short periods of time cause bleeding or if the bleeding persists after removal of the device, stop using it and consult your clinician.

Rapid shrinking of the stoma

If you are wearing the Provox LaryTube because of a rapidly shrinking tracheostoma, discuss with your clinician the possibility of having two devices (of the same size). This would allow you to remove the Provox LaryTube from your stoma and replace it immediately with a clean one. Cleaning of the removed device can then be done after you have inserted the clean LaryTube.

Stoma enlarges during use

If your tracheostoma enlarges after using the LaryTube, ask your clinician if it is possible for you to use two devices (one with a smaller diameter and one with a larger diameter). This allows you to use the larger size LaryTube when the tracheostoma has enlarged, and the smaller size LaryTube to allow the tracheostoma to shrink. It may be necessary to discontinue use of the Tube during the night, so that the tracheostoma can shrink. In that case, discuss alternative options to attach the HME Cassette with your clinician so that you can continue HME Cassette use during the night.

4. Additional information

4.1 Compatibility with MRI Examination

MR-Safe: This device does not contain any metallic elements and has no potential for interaction with the MRI field.

4.2 Ordering information

See end of this Instructions for Use.

4.3 User assistance information

For additional help or information please see back cover of this Instruction for Use for contact information.

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Provox HME System	REF	Rx
Provox XtraMoist HME (30 pcs)	7290	-
Provox XtraFlow HME (30 pcs)	7291	-
Provox FreeHands HME Cassette (20/30 pcs)	7712 / 7711	Rx
Provox FreeHands FlexiVoice Set Plus	7757	Rx
Provox FreeHands FlexiVoice Set	7760	Rx
Provox FreeHands FlexiVoice Light	8161	Rx
Provox FreeHands FlexiVoice Medium	8162	Rx
Provox FreeHands FlexiVoice Strong	8163	Rx

Provox HME System	REF	Rx
Provox FreeHands HME Moist (30 pcs)	8220	Rx
Provox FreeHands HME Flow (30 pcs)	8221	Rx
Provox Micron HME (5/30 pcs)	7247 / 7248	Rx
Provox StabiliBase (15 pcs)	7289	
Provox StabiliBase OptiDerm (15 pcs)	7318	
Provox XtraBase Adhesive (20 pcs)	7265	
Provox Adhesive FlexiDerm Round / Oval / Plus (20 pcs)	7253 / 7254 / 7331	
Provox Adhesive OptiDerm Round / Oval / Plus (20 pcs)	7256 / 7255 / 7332	
Provox Adhesive Regular Round / Oval / Plus (20 pcs)	7251 / 7252 / 7330	

Provox Accessories	REF	Rx
Provox LaryClip	7669	-
Provox TubeHolder	7668	-
Provox TubeBrush	7660-7661	-
Provox Fenestration Punch	7654	Rx
Provox Adhesive Remover (50 pcs)	8012	-
Provox Cleaning Towel (200 pcs)	7244	-
Provox Skin Barrier (50 pcs)	8011	-
Provox Silicone Glue 40 ml	7720	-
Skin Tac Adhesive Barrier Wipes	MS407W	-
Provox ShowerAid	7260	-
Provox BasePlate Adaptor	7263	-
Provox Brush	7204	-

Caution: United States Federal law restricts the sale, distribution and use by, or by order of a physician or a licensed practitioner of devices identified herein with the symbol “Rx”. The prescription-free availability of these products outside the United States may vary from country to country.

This page is deliberately left empty

No. 90734, 2019-06-26

Atos
atosmedical.com

Atos Medical AB

Kraftgatan 8, P.O. Box 183, SE-242 22 Hörby, Sweden

Tel: +46 (0)415 198 00 • info@atosmedical.com